



7975 North Central Drive
Lewis Center, OH 43035
Phone: 740-548-7298 or 800-444-1508
Fax: 740-548-5314 or 800-255-5314
Web: www.globalgilson.com

PRESS RELEASE

Gilson Company, Inc. Announces the New SP-230 Spinning Riffler for Sampling Accuracy and Reliability of the Highest Order [PTXi Expo Booth #2741](#)

Gilson's new SP-230 Spinning Riffler provides the most accurate sampling method available for powders and granular material. Operator error and bias associated with other types of sample dividing are eliminated. The unique design of the SP-230 allows the use of one to sixteen 60ml glass test tubes or 20ml, 40ml, and 67ml vials giving the ability to divide 50-1000ml without an adapter. Unlike most other spinning rifflers the sample vessels are contained within a collection drum during operation to minimize exposure of moving parts and to use as a bulk sample collector. The feeder vibration amplitude and sample collection drum rotation speed are variable and are controlled on the custom touch-screen. The overall design of the unit allows for easy clean-up and maintenance.

For more information, contact Kelly Nardella at 800-444-1508, ext. 153 or knardella@gilsonco.com.

